

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 106

Pellet-burettes, 10 ml, lateral stopc., AS, conf. cert.




Batch: 33-17

Serial:

Measuring Nr. 13987

Batch size: 1

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Measured value:	9,9945 ml	5,9884 ml	1,9894 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	08.09.2017
Balance:	Sartorius - LE244S-OCE

