

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 106

Pellet-burettes, 10 ml, lateral stopc., AS, conf. cert.




Batch: 06-19

Serial:

Measuring Nr. 32442

Batch size: 2

| Volume: | 10,000 ml | 6,000 ml | 2,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value: | 10,0107 ml | 6,0013 ml | 2,0019 ml |
| Deviation: | 0,0033 ml | 0,0026 ml | 0,0028 ml |
| Precision (CV): | 0,0327 % | 0,0426 % | 0,1420 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 30.08.2019 |
| Balance: | Sartorius - LE244S-OCE |

