

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 106

Pellet-burettes, 10 ml, lateral stopc., AS, conf. cert.




Batch: 06-19

Serial:

Measuring Nr. 34631

Batch size: 5

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0059 ml	5,9989 ml	2,0063 ml
Deviation:	0,0070 ml	0,0046 ml	0,0015 ml
Precision (CV):	0,0701 %	0,0763 %	0,0746 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	22.11.2019
Balance:	Sartorius - LE244S-OCE

