

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 108

Burettes, 10 ml, lateral needle valve, AS, conf. cert.

Batch: 22-16




Serial:

Measuring Nr. 12282

Batch size:

6

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	9,9981 ml	5,9946 ml	2,0009 ml
Deviation:	0,0088 ml	0,0047 ml	0,0025 ml
Precision (CV):	0,0883 %	0,0784 %	0,1268 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

25.07.2017

Balance:

Sartorius -
LE244S-OCE

DE-M