

Certificate of performance



Article: 3 122 108

Pellet-burettes, 10 ml, lat. needle valve, AS, conf. cert.




Batch: 02-24

Serial:

Measuring Nr.: 65968

Batch size: 19

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0012 ml	5,9945 ml	2,0005 ml
Deviation:	0,0018 ml	0,0053 ml	0,0065 ml
Precision (CV):	0,0176 %	0,0884 %	0,3251 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	07.02.2024
Balance:	Sartorius - MSE224S-0CE -DU

