

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 255

Pellet-burettes, 25 ml, lat. needle valve, AS, conf. cert.




Batch: 20-23

Serial:

Measuring Nr. 63392

Batch size: 50

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0027 ml	14,9961 ml	5,0123 ml
Deviation:	0,0119 ml	0,0083 ml	0,0032 ml
Precision (CV):	0,0474 %	0,0554 %	0,0635 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	19.06.2023
Balance:	Sartorius - LE244S-OCE

DE-M