

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 258

Pellet-burettes, 25 ml, lat. needle valve, AS, conf. cert.




Batch: 30-19

Serial: -

Measuring Nr.: 36210

Batch size: 17

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9895 ml	14,9962 ml	4,9993 ml
Deviation:	0,0053 ml	0,0100 ml	0,0123 ml
Precision (CV):	0,0213 %	0,0667 %	0,2453 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	17.02.2020
Balance:	Sartorius - MSE224S-0CE -DU

DE-M