

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 258

**Pellet-burettes, 25 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 20-23

**Serial:**

**Measuring Nr.** 70326

**Batch size:** 10

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9899 ml	14,9912 ml	5,0227 ml
Deviation:	0,0059 ml	0,0034 ml	0,0034 ml
Precision (CV):	0,0234 %	0,0225 %	0,0686 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 77

**Date:** 15.11.2024

**Balance:** Sartorius - LE244S-OCE

