

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 504

Pellet-burettes, 50 ml, lateral stopc., AS, conf. cert.




Batch: 29-20

Serial:

Measuring Nr. 59830

Batch size: 10

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0117 ml	29,9982 ml	10,0008 ml
Deviation:	0,0186 ml	0,0089 ml	0,0040 ml
Precision (CV):	0,0373 %	0,0297 %	0,0399 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

25.10.2022

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M