

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 505

**Pellet-burettes, 50 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 30-17

**Serial:**

**Measuring Nr.** 18862

**Batch size:** 10

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0025 ml	30,0313 ml	10,0247 ml
Deviation:	0,0092 ml	0,0051 ml	0,0123 ml
Precision (CV):	0,0184 %	0,0172 %	0,1233 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	15.03.2018
Balance:	Sartorius - LE244S-OCE

