

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 122 508

**Pellet-burettes, 50 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 01-18

**Serial:**

**Measuring Nr.:** 29796

**Batch size:** 20

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0231 ml	30,0374 ml	10,0359 ml
Deviation:	0,0128 ml	0,0038 ml	0,0039 ml
Precision (CV):	0,0256 %	0,0126 %	0,0386 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	16.05.2019
Balance:	Sartorius - LE244S-OCE

