

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 508

Pellet-burettes, 50 ml, lat. needle valve, AS, conf. cert.




Batch: 17-24

Serial:

Measuring Nr.: 69861

Batch size: 20

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0008 ml	30,0220 ml	10,0076 ml
Deviation:	0,0182 ml	0,0113 ml	0,0146 ml
Precision (CV):	0,0364 %	0,0376 %	0,1455 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 77

Date: 29.10.2024

Balance: Sartorius -
LE244S-OCE

DE-M