

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 255

Pellet-burettes, 25 ml, amber glass, AS, conf.cert.



Batch: 25-18

Serial:

Measuring Nr. 23264

Batch size: 14

Volume:	25,0000 ml	5,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	24,9799 ml	5,0008 ml
Deviation:	0,0034 ml	0,0038 ml
Precision (CV):	0,0137 %	0,0761 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	10.08.2018
Balance:	Sartorius - LE244S-OCE

DE-M