

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 256

Burettes, 25 ml, amber glass, AS, conf. cert.

Batch: 22-15




Serial:

Measuring Nr. 12822

Batch size:

1

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Measured value:	25,0268 ml	14,9998 ml	5,0008 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

09.08.2017

Balance:

Sartorius -
LE244S-OCE

DE-M