

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 124 256

**Burettes, 25 ml, amber glass, AS, conf. cert.**

**Batch:** 22-15




**Serial:**

**Measuring Nr.:** 12841

**Batch size:**

1

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Measured value:	25,0289 ml	15,0143 ml	5,0082 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

09.08.2017

**Balance:**

Sartorius -  
LE244S-OCE

