

Certificate of performance



Article: 3 124 258

Pellet-burettes, 25 ml, amber glass, AS, conf. cert.




Batch: 22-13

Serial:

Measuring Nr.: 21422

Batch size: 8

Volume:	25,000 ml	15,000 ml	5,000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9810 ml	14,9954 ml	4,9954 ml
Deviation:	0,0022 ml	0,0040 ml	0,0051 ml
Precision (CV):	0,0087 %	0,0263 %	0,1020 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	04.06.2018
Balance:	Sartorius - LE244S-OCE

