

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 508

Burettes, 50 ml, amber glass, AS, conf. cert.




Batch: 09-16

Serial:

Measuring Nr. 12897

Batch size: 15

| Volume: | 50,000 ml | 30,000 ml | 10,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 49,9719 ml | 29,9896 ml | 10,0069 ml |
| Deviation: | 0,0009 ml | 0,0114 ml | 0,0139 ml |
| Precision (CV): | 0,0017 % | 0,0379 % | 0,1394 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 10.08.2017 |
| Balance: | Sartorius - LE244S-OCE |

