

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 508

Pellet-burettes, 50 ml, amber glass, AS, conf. cert.

Batch: 05-14




Serial:

Measuring Nr. 26687

Batch size:

7

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9888 ml	30,0050 ml	10,0370 ml
Deviation:	0,0138 ml	0,0104 ml	0,0019 ml
Precision (CV):	0,0276 %	0,0347 %	0,0190 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

08.01.2019

Balance:

Sartorius -
LE244S-OCE

DE-M