

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 152 104

**Universal autom. burette, 10 ml, AS, conf. cert.**




**Batch:** 43-19

**Serial:**

**Measuring Nr.** 48520

**Batch size:** 4

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0081 ml	5,9978 ml	1,9962 ml
Deviation:	0,0009 ml	0,0020 ml	0,0035 ml
Precision (CV):	0,0092 %	0,0326 %	0,1728 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	24.08.2021
Balance:	Sartorius - MSE224S-0CE -DU

