

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 152 254

Universal automatic burette, 25 ml, AS, conf. cert.

Batch: 29-20




Serial:

Measuring Nr. 48521

Batch size:

6

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0069 ml	14,9997 ml	5,0039 ml
Deviation:	0,0107 ml	0,0127 ml	0,0106 ml
Precision (CV):	0,0428 %	0,0850 %	0,2117 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

25.08.2021

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M