

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 152 508 FCH

Universal autom. burette, 50 ml, FCH, AS, conf. cert.




Batch: 05-16

Serial: -

Measuring Nr.: 27992

Batch size: 3

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9894 ml	30,0060 ml	10,0265 ml
Deviation:	0,0008 ml	0,0157 ml	0,0229 ml
Precision (CV):	0,0016 %	0,0524 %	0,2288 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	26.02.2019
Balance:	Sartorius - LE244S-OCE

DE-M