

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 012

Mikro burette, 1 ml, str. stopc., glass, AS, conf. cert.

Batch: 13-18




Serial:

Measuring Nr. 35530

Batch size:

2

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0024 ml	0,5981 ml	0,1997 ml
Deviation:	0,0004 ml	0,0010 ml	0,0003 ml
Precision (CV):	0,0424 %	0,1653 %	0,1425 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

15.01.2020

Balance:

Sartorius -
LE244S-OCE

DE-M