

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 022

Micro burette, 2 ml, str. stopc., glass, AS, conf. cert.

Batch: 41-20




Serial:

Measuring Nr. 61471

Batch size:

5

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0020 ml	1,1993 ml	0,4000 ml
Deviation:	0,0009 ml	0,0006 ml	0,0017 ml
Precision (CV):	0,0455 %	0,0537 %	0,4129 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	01.02.2023
Balance:	Sartorius - LE244S-OCE

DE-M