

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 372 022

**Micro burette, 2 ml, str. stopc., glass, AS, conf. cert.**




**Batch:** 38-24

**Serial:**

**Measuring Nr.** 70420

**Batch size:** 10

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
<b>Tolerance:</b>	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
<b>Average value:</b>	1,9949 ml	1,1956 ml	0,4010 ml
<b>Deviation:</b>	0,0024 ml	0,0024 ml	0,0015 ml
<b>Precision (CV):</b>	0,1188 %	0,2024 %	0,3699 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 77

**Date:** 20.11.2024

**Balance:** Sartorius -  
LE244S-0CE

