

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 372 053

**Micro burette, 5 ml, str. stopc., PTFE, AS, conf. cert.**




**Batch:** 10-20

**Serial:**

**Measuring Nr.:** 43211

**Batch size:** 10

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
<b>Average value:</b>	4,9968 ml	2,9937 ml	0,9960 ml
<b>Deviation:</b>	0,0026 ml	0,0011 ml	0,0006 ml
<b>Precision (CV):</b>	0,0529 %	0,0351 %	0,0565 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	77
<b>Date:</b>	11.12.2020
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

