

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 372 102

Micro burette, 10 ml, str. Stopc., glass, AS, conf. cert.




Batch: 46-17

Serial:

Measuring Nr. 24644

Batch size: 1

| Volume: | 10,000 ml | 6,000 ml | 2,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Measured value: | 10,0034 ml | 6,0055 ml | 2,0003 ml |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 4 |
| Date: | 08.10.2018 |
| Balance: | Sartorius - LE244S-OCE |

