

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 012

Micro bur., Bang, 1 ml, lat. stopcock, AS, conf. cert.




Batch: 39-14

Serial:

Measuring Nr. 12109

Batch size: 1

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Measured value:	0,9960 ml	0,5961 ml	0,1979 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	18.07.2017
Balance:	Sartorius - LE244S-OCE

DE-M