

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 383 022

**Micro bur., Bang, 2 ml, lat. stopcock, AS, conf. cert.**

**Batch:** 41-20




**Serial:**

**Measuring Nr.** 51322

**Batch size:** 6

6

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0009 ml	1,2008 ml	0,4015 ml
Deviation:	0,0021 ml	0,0026 ml	0,0012 ml
Precision (CV):	0,1031 %	0,2173 %	0,2994 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

13.01.2022

**Balance:**

Sartorius -  
LE244S-OCE

DE-M