

Certificate of performance



Article: 3 383 025

Micro burettes, 2 ml, lat. needle valve, AS, conf. cert.

Batch: 09-16




Serial:

Measuring Nr. 14108

Batch size:

6

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0021 ml	1,2001 ml	0,3972 ml
Deviation:	0,0013 ml	0,0034 ml	0,0020 ml
Precision (CV):	0,0647 %	0,2815 %	0,5060 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	13.09.2017
Balance:	Sartorius - LE244S-OCE

