

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 025

Micro burettes, 2 ml, lat. needle valve, AS, conf. cert.




Batch: 16-13

Serial:

Measuring Nr.: 24607

Batch size: 6

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0000 ml	1,1998 ml	0,3997 ml
Deviation:	0,0025 ml	0,0018 ml	0,0021 ml
Precision (CV):	0,1254 %	0,1501 %	0,5219 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	08.10.2018
Balance:	Sartorius - LE244S-OCE

