

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 383 052

**Micro burette Bang, 5 ml, lat. stopc., AS, conf.cert.**




**Batch:** 12-20

**Serial:**

**Measuring Nr.:** 46774

**Batch size:** 4

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0049 ml	3,0018 ml	1,0016 ml
Deviation:	0,0008 ml	0,0006 ml	0,0021 ml
Precision (CV):	0,0164 %	0,0217 %	0,2112 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	17.05.2021
Balance:	Sartorius - MSE224S-0CE -DU

