

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 052

Micro burette Bang, 5 ml, lat. stopc., AS, conf.cert.




Batch: 35-21

Serial:

Measuring Nr. 56395

Batch size: 10

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0027 ml	2,9982 ml	1,0016 ml
Deviation:	0,0019 ml	0,0030 ml	0,0020 ml
Precision (CV):	0,0387 %	0,1013 %	0,2011 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	31.05.2022
Balance:	Sartorius - MSE224S-0CE -DU

