

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 055

Micro burettes, 5 ml, lat. needle valve, AS, conf. cert.




Batch: 24-16

Serial:

Measuring Nr.: 25566

Batch size: 10

| Volume: | 5,0000 ml | 3,0000 ml | 1,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value: | 4,9986 ml | 2,9957 ml | 0,9987 ml |
| Deviation: | 0,0038 ml | 0,0017 ml | 0,0010 ml |
| Precision (CV): | 0,0751 % | 0,0577 % | 0,0973 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 09.11.2018 |
| Balance: | Sartorius - LE244S-OCE |

