

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 383 055

**Micro burettes, 5 ml, lat. needle valve, AS, conf. cert.**




**Batch:** 12-20

**Serial:**

**Measuring Nr.** 43329

**Batch size:** 5

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0012 ml	3,0071 ml	1,0020 ml
Deviation:	0,0009 ml	0,0012 ml	0,0003 ml
Precision (CV):	0,0173 %	0,0386 %	0,0304 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	16.12.2020
Balance:	Sartorius - MSE224S-0CE -DU

DE-M