

Certificate of performance



ISO 9001
CERTIFIED

Article: 4 712 322

Pellet-burettes, 25 ml, lat. needle valve, AS, conf. cert.




Batch: 38-15

Serial:

Measuring Nr.: 19026

Batch size: 2

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9926 ml	14,9987 ml	5,0070 ml
Deviation:	0,0064 ml	0,0138 ml	0,0082 ml
Precision (CV):	0,0255 %	0,0918 %	0,1646 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 4

Date: 20.03.2018

Balance: Sartorius -
LE244S-OCE

DE-M