

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 4711632

**Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.**




**Batch:** 50-17

**Serial:** -

**Measuring Nr.:** 24420

**Batch size:** 39

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	49,9879 ml	29,9978 ml	10,0016 ml
Deviation:	0,0110 ml	0,0136 ml	0,0139 ml
Precision (CV):	0,0219 %	0,0454 %	0,1391 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

01.10.2018

**Balance:**

Sartorius -  
LE244S-OCE

DE-M