

Certificate of performance



ISO 9001
CERTIFIED

Article: 4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.




Batch: 32-18

Serial:

Measuring Nr. 33403

Batch size: 50

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0119 ml	29,9990 ml	9,9943 ml
Deviation:	0,0091 ml	0,0129 ml	0,0144 ml
Precision (CV):	0,0182 %	0,0429 %	0,1440 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	26.09.2019
Balance:	Sartorius - LE244S-OCE

DE-M