

Certificate of performance



ISO 9001
CERTIFIED

Article: 4711632

Burettes, 50 ml, lateral valve, Schellb., AS, conf. cert.




Batch: 31-22

Serial:

Measuring Nr. 60725

Batch size: 50

| Volume: | 50,0000 ml | 30,0000 ml | 10,0000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 50,0028 ml | 30,0134 ml | 10,0038 ml |
| Deviation: | 0,0110 ml | 0,0143 ml | 0,0094 ml |
| Precision (CV): | 0,0220 % | 0,0477 % | 0,0944 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|-----------------------------------|
| Operator: | 4 |
| Date: | 13.12.2022 |
| Balance: | Sartorius - MSE224S-0CE -DU |

