

Certificate of performance



Article: 3 715 000 BL

Grad. cylinder, 5:0,1 ml, blue grad., A, conf. cert.




Batch: 50-19

Serial:

Measuring Nr. 35300

Batch size: 120

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	1,0069 ml	3,0059 ml	5,0007 ml
Deviation:	0,0093 ml	0,0109 ml	0,0138 ml
Precision (CV):	0,9290 %	0,3649 %	0,2752 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	19.12.2019
Balance:	Sartorius - LE244S-OCE

