

Certificate of performance



Article: 3 715 000 BL

Grad. cylinder, 5:0,1 ml, blue grad., A, conf. cert.




Batch: 48-24

Serial:

Measuring Nr. 70704

Batch size: 280

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	0,9776 ml	2,9899 ml	5,0012 ml
Deviation:	0,0124 ml	0,0096 ml	0,0097 ml
Precision (CV):	1,2423 %	0,3200 %	0,1933 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	27.11.2024
Balance:	Sartorius - MSE224S-0CE -DU

