

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 000 USP

**Graduated cylinder, 5:0,1 ml, USP, A, conf. cert.**




**Batch:** 34-17

**Serial:**

**Measuring Nr.** 13420

**Batch size:** 90

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	1,0114 ml	3,0070 ml	4,9815 ml
Deviation:	0,0184 ml	0,0162 ml	0,0137 ml
Precision (CV):	1,8438 %	0,5411 %	0,2745 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	24.08.2017
Balance:	Sartorius - LE244S-OCE

