

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 000 USP

Graduated cylinder, 5:0,1 ml, USP, A, conf. cert.




Batch: 04-21

Serial:

Measuring Nr.: 44354

Batch size: 150

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	0,9944 ml	3,0142 ml	5,0135 ml
Deviation:	0,0112 ml	0,0057 ml	0,0092 ml
Precision (CV):	1,1232 %	0,1900 %	0,1830 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	11.02.2021
Balance:	Sartorius - LE244S-OCE

