

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 000

Graduated cylinder, 5:0,1 ml, A, conf. cert.




Batch: 15-18

Serial:

Measuring Nr.: 19865

Batch size: 269

| Volume: | 1,0000 ml | 3,0000 ml | 5,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 0,9823 ml | 2,9857 ml | 4,9865 ml |
| Deviation: | 0,0131 ml | 0,0117 ml | 0,0126 ml |
| Precision (CV): | 1,3070 % | 0,3901 % | 0,2516 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 44 |
| Date: | 12.04.2018 |
| Balance: | Sartorius - LE244S-OCE |

