

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 000

**Graduated cylinder, 5:0,1 ml, A, conf. cert.**




**Batch:** 02-19

**Serial:**

**Measuring Nr.** 26812

**Batch size:** 500

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	1,0072 ml	3,0162 ml	5,0075 ml
Deviation:	0,0163 ml	0,0132 ml	0,0111 ml
Precision (CV):	1,6309 %	0,4391 %	0,2211 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	16.01.2019
Balance:	Sartorius - LE244S-OCE

