

# Certificate of performance



**Article:** 3 715 000

**Graduated cylinder, 5:0,1 ml, A, conf. cert.**




**Batch:** 30-20

**Serial:**

**Measuring Nr.** 39312

**Batch size:** 475

Volume:	1,000 ml	3,000 ml	5,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	1,0027 ml	3,0049 ml	4,9930 ml
Deviation:	0,0202 ml	0,0180 ml	0,0167 ml
Precision (CV):	2,0227 %	0,5991 %	0,3340 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	29.07.2020
Balance:	Sartorius - LE244S-OCE

