

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 000

Graduated cylinder, 5:0,1 ml, A, conf. cert.




Batch: 31-20

Serial:

Measuring Nr.: 39497

Batch size: 100

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	1,0188 ml	3,0077 ml	4,9911 ml
Deviation:	0,0128 ml	0,0166 ml	0,0162 ml
Precision (CV):	1,2773 %	0,5543 %	0,3230 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	07.08.2020
Balance:	Sartorius - MSE224S-0CE -DU

