

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 000

**Graduated cylinder, 5:0,1 ml, A, conf. cert.**




**Batch:** 16-22

**Serial:**

**Measuring Nr.** 54471

**Batch size:** 1200

Volume:	1,0000 ml	3,0000 ml	5,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	0,9903 ml	3,0015 ml	5,0076 ml
<b>Deviation:</b>	0,0120 ml	0,0103 ml	0,0090 ml
<b>Precision (CV):</b>	1,1974 %	0,3447 %	0,1799 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	19.04.2022
<b>Balance:</b>	Sartorius - LE244S-OCE

