

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001 BL

Grad. cylinders, 10:0,2 ml, blue grad, A, conf. cert.




Batch: 06-18

Serial:

Measuring Nr. 18102

Batch size: 500

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9846 ml	5,9966 ml	10,0139 ml
Deviation:	0,0197 ml	0,0169 ml	0,0152 ml
Precision (CV):	0,9860 %	0,2815 %	0,1519 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	20.02.2018
Balance:	Sartorius - LE244S-OCE

DE-M