

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001 USP

Grad. cylinder, 10:0,2 ml, USP, A, conf. cert.




Batch: 04-19

Serial:

Measuring Nr. 27644

Batch size: 100

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	2,0075 ml	6,0026 ml	10,0062 ml
Deviation:	0,0172 ml	0,0148 ml	0,0180 ml
Precision (CV):	0,8606 %	0,2472 %	0,1797 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	12.02.2019
Balance:	Sartorius - LE244S-OCE

