

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 001 USP

**Grad. cylinder, 10:0,2 ml, USP, A, conf. cert.**




**Batch:** 46-20

**Serial:**

**Measuring Nr.:** 42799

**Batch size:** 280

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9998 ml	6,0059 ml	10,0105 ml
Deviation:	0,0147 ml	0,0134 ml	0,0159 ml
Precision (CV):	0,7330 %	0,2230 %	0,1588 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	25.11.2020
Balance:	Sartorius - MSE224S-0CE -DU

